

Revised August 2009



HONORS LAB 13b: Iodine Clock


Aim To determine the orders of reaction in a specific chemical reaction

Apparatus Graduated cylinders (various sizes), 400 mL beakers, stopwatch

Chemicals 0.050 M KI, 0.050 M Na₂S₂O₃, starch solution, buffer solution (pH = 4.7), 0.80 M H₂O₂, deionized water

Introduction

The chemical reaction to be investigated is given below



Method

1. Obtain three, 400 mL beakers and label them A-C.
2. To the relevant beaker, add deionized water and the KI solution in the quantities specified in the table below.


Beaker	Volume of deionized water in mL	Volume of 0.05 M KI in mL
A	125	25.0
B	100.	50.0
C	115	25.0

3. To each beaker, also add 5.0 mL of the starch solution to act as an indicator, 30.0 mL of the buffer solution and 5.0 mL of the Na₂S₂O₃ solution.
4. Starting with beaker A, quickly add 10.0 mL of the H₂O₂ solution while simultaneously starting the stopwatch. Stop the stopwatch when the intense blue/black color is observed. Record the time taken for the blue/black color to appear in the results table.

Beaker	Volume of 0.80 M H ₂ O ₂ in mL
A	10.0
B	10.0
C	20.0

5. Repeat the procedure in # 4 for beakers B and C on each occasion using the volume of H₂O₂ solution specified in the table above.

Revised August 2009



Results

Given that the TOTAL volume of solution in each beaker is 200. mL, fill in the table below.

Beaker	[I ⁻]	[H ₂ O ₂]	Time taken in s
A			
B			
C			

